

Amendments to the Specification:

Please replace paragraph [0015] on page 4 with the following amended paragraph:

A balloon catheter capable of both fast and simple guidewire and catheter exchange is particularly advantageous. A catheter designed to address this need is sold by Medtronic Vascular, Inc. of Santa Rosa, Calif. under the trademarks MULTI-EXCHANGE, ZIPPER MX, ZIPPER, MX and/or MXII (hereinafter referred to as the "MX catheter"). An MX catheter is disclosed in U.S. Pat. No. 4,988,356 to Crittenden et al.; ~~eo-pending U.S. patent application Ser. No. 10/116,234, filed Apr. 4, 2002 U.S. Patent No. 6,800,065 to Duane et al.;~~ ~~eo-pending U.S. patent application Ser. No. 10/251,578, filed Sep. 18, 2002 U.S. Patent Appl. Publ. No. 2004/0059369 A1 to Duffy et al.;~~ ~~eo-pending U.S. patent application Ser. No. 10/251,477, filed Sep. 20, 2002 U.S. Patent No. 6,905,477 to McDonnell et al.;~~ ~~eo-pending U.S. patent application Ser. No. 10/722,191, filed Nov. 24, 2003 U.S. Patent Appl. Publ. No. 2004/0260329 A1 to Gribbons et al.;~~ and ~~eo-pending U.S. patent application Ser. No. 10/720,535, filed Nov. 24, 2003 U.S. Patent No. 6,893,417 to Gribbons et al.~~, all of which are incorporated by reference in their entirety herein.

Please replace paragraph [0018] on page 5 with the following amended paragraph:

A clinician may wish to perform fast and simple guidewire and catheter exchanges while maintaining a guidewire ~~guidewire~~ fully within a catheter as in a conventional OTW catheter. An alternative form of guide member that allows that capability (hereinafter referred to as the "grabber") is disclosed in ~~eo-pending U.S. patent application Ser. No. 10/226,789, filed Aug. 21, 2002 U.S. Patent Appl. Publ. No. 2004/0039372 to Carmody et al.~~, that is incorporated by reference in its entirety herein. The grabber is similar to the guide member described above in that it is slidably coupled to a MX catheter shaft. However, the grabber does not allow a guidewire to enter or exit the MX catheter anywhere along the length of the catheter shaft. Instead, the grabber allows a clinician to apply a clamping force on a guidewire within the catheter shaft allowing him to directly manipulate the position of the guidewire within the catheter shaft.